




Echem Analyst™ Software Operator's Guide



Copyright © 2019 Gamry Instruments, Inc.
Revision 1.2
April 1, 2019
Gamry Instruments P/N 988-00016

If You Have Problems

Please visit our service and support page at www.gamry.com/service-support/. This page contains information on installation, software updates, and training. It also contains links to the latest available documentation. If you are unable to locate the information you need from our website, you can contact us via email using the link provided on our website. Alternatively, you can contact us one of the following ways:

Internet	www.gamry.com/service-support/
Telephone	(215) 682-9330 9:00 AM–5:00 PM, US Eastern Standard Time (877) 367-4267, Toll Free US and Canada Only

Please have your instrument model and serial numbers available, as well as any applicable software and firmware revisions.

If you have problems in installation or use of a system containing Echem Analyst software, please try to call from a phone next to your computer, where you can type and read the screen while talking to us.

We will be happy to provide a reasonable level of free support for registered users of the Echem Analyst software. Reasonable support includes telephone assistance covering the normal installation, use, and simple customization of a computerized system containing Echem Analyst software connected to a Windows®-compatible computer.

A service contract that extends the software update period is available at an additional charge. Software updates **do not** include software enhancements offered to our customers at additional cost.

Enhancements to the Echem Analyst software that require significant engineering time on our part can be performed on a contract basis. Contact us with your requirements.

Copyrights

Echem Analyst™ Software Operator's Guide copyright © 2008–2019, Gamry Instruments, Inc., all rights reserved.

Framework copyright © 1989–2019, Gamry Instruments, Inc., all rights reserved.

Framework™ and Gamry™ are trademarks of Gamry Instruments, Inc.

License Agreement

Gamry Instruments, Inc. grants you, the Purchaser, a non-exclusive license to use the software programs included in this package, including updates which may later be provided by Gamry Instruments at additional cost. This license is subject to the terms and restrictions set forth in this Agreement.

License

The software programs and the Operator's Manuals are protected by both United States copyright law and international treaty provisions. Accordingly, these software programs remain the property of Gamry Instruments, Inc. as its sole property, subject to the rights granted to you in this License Agreement.

You may

- Make archival copies of the software programs for backup purposes only. All backup copies must include the Gamry Instruments, Inc. copyright notice.
- Transfer the programs and Gamry potentiostat from computer to computer within your facility.
- Operate the data-acquisition (Framework™) and data-analysis portion (Echem Analyst™) of the programs on an unlimited number of computers within your facility.
- Transfer the programs, with a copy of this license, to another party, provided the other party agrees to accept the terms and conditions of this Agreement. If you transfer the programs, you must at the same time transfer all copies of the programs to the same party or destroy any copies not transferred. If you transfer possession of any copy of the software programs, in whole or in part, to another party, your license is automatically terminated.

You may not

- Sublicense, rent or lease the software programs.
- Use this program with a potentiostat not specifically authorized by Gamry Instruments, Inc.
- Modify the program and/or potentiostats to defeat the authorization code security scheme.
- Reverse engineer, decompile, or disassemble any program module or programmable logic device.
- Copy, use, modify, or transfer the software programs or Operator's Manuals or any copy, in whole or in part, except as expressly provided in this Agreement.

Term

This license is effective until terminated. You may terminate it by destroying or returning the programs and Operator's Manuals and any copies thereof. This license will terminate if you fail to comply with any term or condition of this Agreement.

Governing Law

1. This License Agreement shall be governed by the laws of the Commonwealth of Pennsylvania. If any provision of this License Agreement is found void, invalid, or unenforceable, it will not affect the other terms and conditions of this Agreement, which remain valid and enforceable.

Third-Party / Use Sublicenses

2. Gamry software may incorporate Third-Party or Use Sublicenses (such as Microsoft VBA). In this case, the terms and conditions of their License Agreements (available during Setup) will apply.

Disclaimers

Gamry Instruments, Inc. cannot guarantee that the Echem Analyst software will work with all computer systems, operating systems, and third-party software applications, hardware, or software.

The information in this guide has been carefully checked and is believed to be accurate as of the time of release. However, Gamry Instruments, Inc. assumes no responsibility for errors that might appear.

No part of this document may be copied or reproduced in any form without the prior written consent of Gamry Instruments, Inc.

Windows®, Excel®, and Powerpoint® are registered trademarks of Microsoft Corporation. Origin® is a registered trademark of OriginLab Corporation. DigiSim® is a registered trademark of Bioanalytical Systems, Inc.

Echem Analyst™ Software

Echem Analyst™ Software.....	5
Introduction to This Guide	6
General Information and Overview	7
Installation	7
File Formats	7
To Open a Gamry Data File.....	7
Working with Plots in Echem Analyst.....	9
Changing the Axes on a Plot (the <i>Curve Selector</i>):	12
Selecting Portions of a Curve for Analysis.....	12
Cutting and Pasting Images and Data.....	14
Getting On-Line Help and Version Number:	16
Common Tools and Tabs.....	17
Accessing Common Tools	17
List of <i>Common Tools</i>	17
Experimental Setup.....	19
Experimental Notes.....	20
Hardware Settings.....	21
Open Circuit Voltage (Corrosion Potential) Data.....	23
Analysis of Cyclic Voltammetry Data.....	24
Cyclic Voltammetry Special Tools	24
Integrating the Voltammogram.....	26
Modeling Polarization Resistance Data	28
Polarization Resistance Special Tools	28
Finding the Polarization Resistance	28
Modeling Potentiodynamic (Tafel) Data.....	30
Tafel Fit.....	30
E Log I Fit.....	31
Modeling EIS Data	32
Bode and Nyquist Plot View	32
EIS Special Tools	33
Fitting the Data to the Equivalent-Circuit Model	37
Appendix.....	40
Headings in Data-File Columns.....	40
Current Conventions According to Framework™ and Echem Analyst.....	41
To Edit Visual Basic Scripts:.....	41
Simulating an EIS Curve	42

Introduction to This Guide

Echem Analyst™ software is Gamry Instruments' dedicated data-analysis program, the companion to Gamry Instruments' data-acquisition program called Framework™. Data files generated by experiments in Framework software then can be analyzed in the Echem Analyst. The Echem Analyst is a single program that runs data-analysis for all types of experiments, such as those used in DC Corrosion, EIS, and Physical Electrochemistry.

The Echem Analyst is designed with the specific functions to make data analysis as straightforward as possible. This manual explains the most common analysis routines. The tools discussed here in the examples are common to many data files created by other experiments. This document is a guide, and is not intended to have the same scope as the on-line help or a complete operating manual. In order to create a concise document, we assume you have a working knowledge of Windows®-based applications. We ignore details about common functions, such as opening, saving, and closing files, so as not to obscure the goal of this guide.

This textbox indicates a helpful hint to know about Echem Analyst.



General Information and Overview

Installation

Echem Analyst installs separately from other Gamry software. If Echem Analyst is not installed yet, you can find it on the CD-ROM, or—if you already own one of our potentiostats—on our website at www.gamry.com.

You may install copies of the Echem Analyst on multiple computers. Often users prefer the convenience of performing data-analysis at an office workstation, rather than the laboratory setting.

File Formats

Gamry data files acquired using Framework software have the extension *.DTA. DTA files are ASCII text, and therefore you can open them directly into various programs, such as Excel® or Origin®. When DTA files are opened in Echem Analyst, then saved, their extension becomes *.GData. Gdata (Gamry Data) files include information on curve-fits and graphing options, thus Gdata files are only viewable in Echem Analyst.



Do not delete your DTA files! They are the raw data and may need to be reloaded for certain analyses, such as area normalization.

To Open a Gamry Data File

There are several different methods to open data files in the Echem Analyst:

3. Launch the Echem Analyst icon on your desktop. Then use the **File** > **Open** function.
4. Use the link on your desktop to open the **My Gamry Data** folder. Double-click on the data file. You may have to instruct your computer to associate the *.DTA extension with the Echem Analyst program.
5. There are two quick ways to open a recent Gamry Data File.
 - a. A recently generated file can be opened using the hotlink in the *Analysis* menu in the Gamry Framework. (The last eight generated data files are listed there for quick access.) The Echem Analyst automatically launches and opens your selected data file.
 - b. A recently opened file in the Echem Analyst is shown at the bottom of the **File** menu. This is similar to how other Windows®-based programs display links to *Most Recently Used* documents.




By default, files acquired in the Framework are saved into the **My Gamry Data** folder. A shortcut for **My Gamry Data** installs on the Windows® desktop. You can change this default under **Tools > Options**, which opens the **Gamry Analysis Framework Options** window. Choose the **General** tab, and change the **Path** for each type of data file as desired.



Don't change the directory for *Analysis Script* Files. These are the VBA programs that do the actual analysis.

The data set appears in the main window. The menu items, tabs, and toolbar are adjusted for the particular type of data set you chose. In the example below, a Potentiostatic EIS data set is shown:



Note the tab-based display. The **Experimental Setup** tab displays all the information from the parameters used to run the experiment, such as voltage, time, etc. The **Experimental Notes** tab stores any notes written into the setup screen in Framework. The **Open Circuit Voltage** tab shows the voltage measured during the initial delay of the experiment. The **Hardware Settings** tab records information on the filters, ranges, gains. Additional information on date of last calibration, software version, etc. is also stored here.


Working with Plots in Echem Analyst

Introduction

Echem Analyst boasts a number of graphical tools to help you get the most information out of your data. After you open a data set, these tools appear in the toolbars immediately above the plot:

In the data file, Framework writes a line that indicates the type of experiment used to generate that file. Echem Analyst displays both general and specific menus containing the analysis routines pertinent to your experiment.

Toolbars:



The main toolbars are:


General
Toolbar

General functions for replotting and printing in various formats


















Selection
Toolbar

Tools to select and view data points













The following charts are references for buttons on the default toolbars. Descriptions of the most commonly used functions are highlighted in blue.

General Toolbar *Functions*

Button	Name	Action
	Copy to clipboard	Copy the selection to the Windows® clipboard. Can paste directly in Microsoft programs for reports or presentations.
	Gallery	Choose, via the dropdown menu, from scatter (no line), line, curve, and steps between data points
	Color	Choose the color of the selection from the dropdown menu. To change the color menu, use the <i>Palette</i> button on the PaletteBar .
	Vertical Grid	Toggle between showing and hiding vertical grid lines on the plot
	Horizontal Grid	Toggle between showing and hiding horizontal grid lines on the plot
	Legend Bar	Toggle between showing and hiding a legend bar underneath the plot
	Data Viewer	Toggle between showing and hiding numerical data to the left of the plot
	Properties...	Open the GamryChart Properties window, so that you can adjust effects, colors, markers, 3-D effects, lines, etc.
	3D/2D	Toggle between two-dimensional and three-dimensional graphing
	Rotate	Rotate the three-dimensional graph. Only active if the graph is 3D.
	Z-clustered	Offset two data sets so that they can be distinguished within one plot. Only operates in 3D mode.
	Zoom	Zoom in on a selected region. Also open a zoom slider at the bottom of the graph for continuous adjustment of zoom.
	Print preview	Open the Page setup window to adjust orientation of plot and printer margins
	Print	Print the plot
	Tools	Open a dropdown menu, for choices of various toolbars and viewers to appear on the screen

Selection Toolbar *Functions*

Button	Name	Action
	Show curve selector	Open the Curve Selector area to the right of the plot, so that you can choose which data are used as the x-, y-, or y2-coordinate, and which curve is the active trace.
	Select x region	Select a desired region of the plot across the x-axis. Commonly used to specify a region for Quick-Integrate .
	Select y region	Select a desired region of the plot across the y-axis. Commonly used to specify a region for Quick-Integrate .
	Select Portion of Curve using the Mouse	Left-click on the active trace using the mouse to select a section of the curve
	Select Portion of Curve using the Keyboard	Open an area to the right of the plot, in which you can choose a segment of the trace numerically. See below for more details.
	Draw Freehand Line	Draw a line on the plot
	Mark Found Peaks	Place a tag on peaks that the software finds. A portion of the curve must be selected first.
	Apply Template	Open the Apply User-Defined Chart Template window, and choose a previously created template to apply to the plot
	Save Template	Open the Save User-Defined Chart Template window, and save the template
	Show Disabled Points	Show data points not being used in the plot

Changing the Axes on a Plot (the *Curve Selector*):

To choose a different variable plotted on an axis, use the **Curve Selector** button




as follows:

(The example shown below is a Differential Pulse Voltammetry plot.)

1. With the plot open and displayed on the screen, click the **Curve Selector** but-

ton



The **Curve Selector** area appears on the right side of the window.

2. Choose which trace is active by clicking the drop-down menu in the **Active Trace** area.

The Active Trace is the data series on which the analysis will be performed. Use this if multiple files or cycles are displayed on the graph.

3. Choose which trace is visible on the plot by activating the checkbox next to the desired trace(s) in the **Visible Traces** area.

Visible Traces also contains any data fits that are performed.


4. Choose which variable is plotted on the x-axis by highlighting the variable in the **X-Axis** column.
5. Choose which variable is plotted on the y-axis by highlighting the variable in the **Y-Axis** column.
6. Choose which variable is plotted on the second y-axis by highlighting the variable in the **Y2-Axis** column.




If there is a data column graphed on the **Y2-Axis**, those data appear in a different color and a different scale.

Selecting Portions of a Curve for Analysis

For certain types of analysis, you must select a region of the curve, for example, within the **Peak Find** function in Cyclic Voltammetry or **Tafel Fit** function in Potentiodynamic. You can select regions by mouse or keyboard.

1. Left-click the mouse on the **Mouse** button  in the **Selection toolbar**.
2. Use the left mouse-button to select each endpoint of the curve. Each endpoint is marked with a blue cross. The selected portion of the curve is shown as a thick blue line. (In the figure below, the color of the data was changed to red for contrast to the selected region).



3. Another click on the **Mouse** button clears the selected region, and readies the graph for a different region to be selected.


Cutting and Pasting Images and Data

Many users want to present, publish, or otherwise share their data and charts from the Echem Analyst.

To create a bitmap image of the graph,


As a Bitmap:

1. In the General Toolbar, click the **Copy to Clipboard** button .
2. In the drop-down menu select **As a Bitmap**.



3. A bitmap image of the graph enters the clipboard. You can paste this bitmap into a presentation program such as Word® or Powerpoint®.


This is a quick and easy way to import a picture of the graph for a presentation or report.



As Text:


Because Gamry Data Files are ASCII text, they can be opened easily in other graphing programs, such as Excel® or Origin®. Right-click on the DTA file, select **Open With...**, and select your desired program. These programs, however, do not contain fitting routines specific to the analysis of electrochemical data. This **As Text** feature lets you fit the data in Echem Analyst and then copy and paste the data and fit into another graphing program.

This is a quick and easy way to import both the data and the fit into another graphing program.



Plotting Conventions

By right-clicking the mouse on a non-zero value on an axis, you can choose to show that axis in **Log Scale** (logarithmic) or linear scale (when **Log Scale** is inactive), or to reverse the direction of the numbers (when **Reversed** is active).



Or you can use the **Transform Axes** selection (if available) under the **Common Tools** menu.



Default plotting of graphs is auto-scale. Therefore, note the y-axis's scale when a plot first appears. If bad data points obscure your data because of auto-scaling, you can choose to disable and hide those offending points.

Getting On-Line Help and Version Number:


In the toolbar, choose **Help**.


- a. Click **Gamry Help** to obtain information about various commands and functions within Echem Analyst.

A separate **Gamry Echem Analyst Help**

window appears. You can find much information

about the details of Echem Analyst here, such as plotting and analysis.





On-line help is a great resource for more involved questions. **Help** is divided up according to software package.

- b. Click **About the Gamry Echem Analyst** to view the software version number.

Common Tools and Tabs

While each type of experimental data has its own method and parameters, there are certain commands that are common to many analyses. This section shows you these **Common Tools**.


Accessing Common Tools

1. Open a dataset.

In the toolbar, the function **Common Tools** appears.

2. Choose **Common Tools**.

A drop-down menu appears.



3. Select the desired command.

In this example, chronoamperometry's **Common Tools** includes three commands: **Add I Constant**, **Linear Fit**, and **Smooth Data**.




The list of **Common Tools** varies depending upon the type of experiment.

List of Common Tools

Command	Type of experiment	Result
Add E Constant	Cyclic Voltammetry, DC Voltammetry, Differential Pulse Voltammetry, Galvanic Corrosion, Normal Pulse Voltammetry, Pitting Scan, Polarization Resistance, Potentiodynamic Scan, Square-Wave Voltammetry	Adds a constant potential to all voltages in the plot. Used to easily convert between different Reference Electrode's scales.
Add I Constant	Chronoamperometry, Chronopotentiometry, Cyclic Voltammetry, Galvanic Corrosion, Pitting Scan, Polarization Resistance, Potentiodynamic Scan	Adds a constant value to all currents in the plot.
C from CPE, omega(max)	Potentiostatic EIS, AC Voltammetry, Mott-Schottky	Calculates capacitance from previously fit CPE values

		and data from the Nyquist plot.
C from CPE, R(parallel)	Potentiostatic EIS, AC Voltammetry, Mott-Schottky	Calculates capacitance from previously fit CPE and fit R data.
Linear Fit	Chronoamperometry, Potentiostatic EIS, AC Voltammetry, Chronocoulometry, Chronopotentiometry, Cyclic Voltammetry, DC Voltammetry, Differential Pulse Voltammetry, EMF Trend, Galvanic Corrosion, Mott-Schottky, Normal Pulse Voltammetry, Polarization Resistance, Potentiodynamic Scan, Square-Wave Voltammetry	When a region of the plot is selected, fits the data to $y = mx + b$.
Post-Run iR Correction	Cyclic Voltammetry, Polarization Resistance, Potentiodynamic Scan	Corrects a previously run scan for voltage-drop caused by iR .
Smooth Data	Chronoamperometry, Chronopotentiometry, Cyclic Voltammetry, DC Voltammetry, Differential Pulse Voltammetry, EMF Trend, Galvanic Corrosion, Normal Pulse Voltammetry, Pitting Scan, Polarization Resistance, Potentiodynamic Scan, Square-Wave Voltammetry	Smooths the data. Useful for locating peaks in regions of high data-density.
Transform Axes	Galvanic Corrosion, Pitting Scan, Polarization Resistance, Potentiodynamic Scan	Changes x- and y-axes from linear to logarithmic, etc.

Experimental Setup




This particular **Experimental Setup** tab is from a Cyclic Voltammetry experiment. This example has many of the same parameters as other experiments. It shows:

Initial E, Scan Limit 1, Scan Limit 2, Final E	The potentials defining the waveform, and whether measured vs. a reference electrode (E_{ref}) or the open circuit potential (E_{oc}).
Test Identifier	Read from the Framework Setup. This field also becomes the default title of the plot.
Time	Time the experiment was started
Scan Rate	How fast (in mV/s) the scan was taken
Step Size	The interval between potentials
Electrode Area	The size of the electrode
Equil. Time	How much time was spent letting the electronics settle before the scan was started
I/E Range Mode	Automatically adjusted or fixed I/E (Current) Range mode.
Max Current	The current value that sets the I/E Range in Fixed Mode and determines the range in which to start in Auto Mode
Conditioning	Whether off or on, for how long, and under what potential. This Potential is vs. Reference.
Init. Delay	Whether off or on. This is when the E_{oc} is measured.
Cycles	Number of how many voltammetry cycles were run
IR Comp	If IR Compensation was used and the mode.
Open Circuit	The value of the open-circuit voltage (Corrosion Potential). It is the value of the last point in the Initial Delay.
Sampling Mode	Data-acquisition mode (for Reference family Potentiostats)

Experimental Notes

Click the **Experimental Notes** tab:

Any notes entered in the Framework software are automatically displayed here. You may enter any additional comments about the experiment in the **Notes...** field.



This is a version of a modern laboratory notebook. Enter as many details about your experiment as you can. Information here can help you avoid having long strings of descriptive file names.

Hardware Settings




This section documents the hardware settings that were used when the experiment was run, e.g., everything from the offsets, filters, and gains to the last time the potentiostat was calibrated.

This information is used primarily by Gamry Instruments' Technical Support staff to help troubleshoot. Gamry Instruments determines default settings based on experience. Advanced users can adjust these settings manually before the experiment is run.

For DC Corrosion experiments, the **Hardware Settings** are set in the experiment code. For Physical Electrochemistry experiments, you have access to these features through the **Advanced Panel**, but we recommend that only advanced users make changes to these settings. Consult **Help** or Gamry Instruments' Technical Support for advice.

Click the **Hardware Settings** tab:



The hardware settings displayed here are:


Potentiostat	Shows the potentiostat's label
Control Mode	How the experiment was controlled
Control Amp Speed	Shows the speed of the control amplifier
I/E AutoRange	Shows if the I/E autorange function was enabled
I _{ch} AutoRange	Shows if the I _{ch} autorange function was enabled
I _{ch} Range	Shows the I _{ch} range (gain). 3 volts = ×1 Gain.
I _{ch} Filter	Shows the I _{ch} cut-off filter frequency
I _{ch} Offset Enable	Shows if I _{ch} Offset was enabled
I _{ch} Offset	Shows the I _{ch} offset voltage
Positive Feedback IR Comp	Shows if the IR positive feedback was enabled
I/E Range Lower Limit	Shows the lowest available I/E Range available to use in this experiment
A _{ch} select	Shows the input connector for A _{ch}
DC Calibration Date	Shows the date of last DC calibration
Framework Version	
Pstat Model	Gives the model number of the potentiostat

Current Convention	Shows which currents are positive
I/E Stability	Shows the I/E stability speed
I/E Range	Shows the I/E (or current) range used
Vch AutoRange	Shows if V_{ch} autoranging is enabled
Vch Range	Shows the maximum value for V_{ch}
Vch Filter	Shows the V_{ch} cut-off filter frequency
Vch Offset Enable	Shows if V_{ch} Offset was enabled
Vch Offset	Shows the I_{ch} offset voltage
Positive Feedback Resistance	Shows the positive feedback resistance applied to the system
Ach Range	Shows the voltage range of the auxiliary channel
Cable ID	(for Reference family potentiostats only.) Gives the type of cable connected to the system
AC Calibration Date	Shows the date of last AC calibration
Instrument Version	Shows the Firmware Version of a Reference family potentiostat

Detailed explanations of these parameters are beyond the scope of this guide.

Open Circuit Voltage (Corrosion Potential) Data


Click the **Open Circuit Voltage** tab:



Because default plotting of graphs is auto-scale, please note the y-axis's scale when the **Open Circuit Voltage** first appears.

Analysis of Cyclic Voltammetry Data


This is a sample cyclic voltammetry file that installs in \My Gamry Data\ when Framework software installs.



Cyclic Voltammetry Special Tools

This menu analyzes the cyclic voltammetry data.

1. In the main menu, choose **Cyclic Voltammetry**.
A drop-down menu appears.
2. Choose the desired tool:



Tool	Function	Notes
Min/Max	Finds the minimum and maximum currents and voltages within the dataset. Results appear in a window below the plot.	
Quick Integrate	Integrates to find the total charge. Results appear in a window below the plot.	For multi-cycle CV experiments
Integrate	Integrates over a specified portion of the plot to find the total charge.	Portion of the curve must be selected
Region Baselines	Defines a line as the baseline for a specified region.	Region must be selected
Clear Regions	Clears all baselines from the dataset.	Region must be selected
Normalize by Scan Rate	Normalizes the dataset based on the scan rate.	
Normalize by Square Root of the Scan Rate	Normalizes the dataset based on the square-root of the scan rate.	
Peak Find	Finds peaks within a specified region of the dataset.	Portion of the curve must be selected
Clear Peaks	Clears all peaks found within the dataset.	Peaks must be identified
Automatic Baseline	Finds the baseline automatically.	Peaks must be identified
Peak Baselines	Defines a line as a baseline for a specified peak.	Peaks must be identified
Clear Lines	Clears all lines from the dataset.	Lines must be associated with graph
Delta Ep	Finds the potential difference between two peaks.	Peaks must be identified
Subtract Background from File	Subtracts a background amount from the dataset.	
Export to DigiElch Options	Exports the file to a DigiElch compatible format. Changes units and grids for plotting the data.	

Integrating the Voltammogram

All integration methods integrate current versus time to get the total charge. There are two different ways to integrate under a curve with Echem Analyst.

Quick Integrate

Quick Integrate breaks the data into “curves.” Each curve is integrated to a zero current. **Quick Integrate** integrates the entire area of each curve, unless an area is specified using the x-region icon.

Integrate

Integrate requires you first to select a portion of the curve. (See how to select a portion of the curve in the “Starting Echem Analyst” chapter.) After an integration is performed, you can change the baseline from the default 0 A to another line, either a line that you draw, or an **Automatic Baseline**.

3. Open the data file.

1. Select the **Draw a Freehand Line** button.



2. Left-click and hold on the graph to place an anchor point. Holding down the mouse button, extend the line with the mouse. Move or extend the line as you wish.

Directions to accept the line are printed at the bottom of the window.

3. Right-click the mouse on the line and either **Accept** or **Delete**.

After you accept the line, it turns from dashed to solid.



4. Select the portion of the curve to integrate. This function is described in detail earlier.



5. Select **Integrate** from the Cyclic Voltammetry menu.


This integrates the section between the curve and the zero A line.

6. To change the baseline to the desired user-drawn line, select **Region Baselines** from the Cyclic Voltammetry menu.




The **Region Baseline Settings** window appears.

7. Select the desired **Baseline** from the available lines. You may draw multiple lines from which to choose.



The integrated region moves from the default 0 A baseline to the user-drawn line.




Modeling Polarization Resistance Data

Polarization Resistance Special Tools

This menu analyzes the polarization resistance data.

1. In the main menu, choose **Polarization Resistance**.

A drop-down menu appears.



2. Choose the desired tool:


Tool	function
Quick Integrate	Integrates to find the total charge. Results appear in a window below the plot.
Min/Max	Finds the minimum and maximum currents and voltages within the dataset. Results appear in a window below the plot.
Polarization Resistance	Within a selected portion of the curve, finds the polarization resistance.
Options	Changes units and grids for plotting the data.

Finding the Polarization Resistance


Method 1: Manual Entry of the Tafel Constants

1. Select the desired portion of the curve. (See section "Selecting Portions of a Curve for Analysis")
2. In the main menu, choose **Polarization Resistance**.

A drop-down menu appears.




3. Choose **Polarization Resistance**.
The **Polarization Resistance** window opens.
4. In the **Seed Values** area, enter anodic (**Beta A**) and cathodic (**Beta C**) Tafel constants.
5. Click the **Calculate** button.
The calculated **Corrosion Rate** appears in a window below the plot.



Method 2: Automatic Selection of Voltage Region

Echem Analyst offers another way to select automatically the voltage region over which this analysis is done.


1. In the Polarization Resistance menu, choose Options.



The **Polarization Resistance Options** window opens.

2. Select the **Automatic** radio button, specify the region around E_{corr} to use, and click the **Save as Defaults** button. You are prompted directly for Tafel constants when a polarization resistance file is opened.

This is how Gamry Instruments' RpEc Trend experiments calculate corrosion rate.




Modeling Potentiodynamic (Tafel) Data

A Tafel experiment is also a very popular electrochemical corrosion technique. The following analysis is performed on the sample Potentiodynamic data file.

Tafel Fit


3. Select the region over which to perform the Tafel fit.
This region must encompass the E_{corr} (open-circuit potential).
4. Select **Tafel Fit** from the **Potentiodynamic** menu:



1. A **Tafel Fit** window appears where you may input seed values optionally for the fit. The better the information we provide the fitting routine, the more likely it will be able to generate an acceptable fit.



If you have reasonable starting parameters for the fit, input them in the **Seed Values** area, and check the **Use Seed Values** checkbox. If you do not have any confidence at all in your range of parameters, do not check the **Use Seed Values** checkbox.

We recommend using the seed values supplied by the Echem Analyst.



2. Click the **Calculate** button.
When you click the **Calculate** button, the changes can be subtle. The following events occur:
 - The parameters in the **Tafel Fit** window become the fit parameters.
 - A fit line is displayed on the graph.

- A new **Tafel** tab is created (to the right of the **Hardware Settings** tab) that holds the information about the fit.




E Log I Fit

The **E Log I** fit is a useful fit if you want to fit the data one branch (anodic or cathodic) at a time. This can be important if one branch doesn't show linear behavior, but the other does.


The fit is called **E Log I** because of the semi-logarithmic nature of a Tafel plot. The x-axis is the logarithm of current, while the y-axis is potential on a linear scale.

Method

5. Select a portion of the curve.
Here you need only the linear section of one of the branches. This selection does not include **Ecorr** [Eoc (open-circuit potential)].
1. In the **E Log I Fit** window, enter an approximate value for **Ecorr**.
2. Click the **Calculate** button.
A single branch of the Tafel data is fit. The fit is shown on the graph, and the results of the fit are contained in a new **E Log I** tab.



You can run a Polarization Resistance fit on these Potentiodynamic data, if the axes of current are changed to the linear scale. Generally we suggest running a separate experiment on a new sample of the same material because of the more-polarizing, more-destructive nature of the Potentiodynamic experiment.




Modeling EIS Data


The data-analysis features shown here are common to many of the AC-based techniques. By far the most popular type of AC experiment is Potentiostatic EIS.

Bode and Nyquist Plot View

Click the **Bode** tab or the **Nyquist** tab of the plot you prefer to work with. All fits are displayed on both the Bode and Nyquist plots. Because they are different representations of the same data, the fit results are identical.



Bode plot



Nyquist plot



EIS Special Tools


EIS data-analysis uses an equivalent-circuit approach. This menu creates and runs fits for EIS data. Commands in this menu allow you to build an equivalent-circuit model in the **Model Editor**, then fit that model to your data. This menu also lets you run advanced procedures, such as **Subtract Impedance**, and run Kramers-Kronig transforms.

6. In the main menu, choose **Impedance**.

A drop-down menu appears.

7. To create or edit an equivalent circuit, choose **Model Editor**.


The **Impedance Model Editor** window appears. See the next page for how to use it.



1. Choose your fitting method.
 - To fit the data using the Levenberg-Marquardt method, choose **Fit A Model (Levenberg-Marquardt Method)**.
 - The **Select Model File** window opens.
 - Choose the appropriate model file, and click the **OK** button.
 - To fit the data using the Simplex method, choose **Fit A Model (Simplex Method)**.
 - Simplex method weighs the user's seed values less. We recommend using the Simplex method.
2. To subtract an impedance from the data, choose **Subtract Impedance...**
The **Impedance Subtraction** window appears.
Choose:
 - Element Choose a circuit element from the drop-down menu.
 - Model Browse for a previously defined model.
 - Spectrum Browse for a data-set.Click the **Close** button.
3. To use the Kramers-Kronig method, choose **Kramers-Kronig**.
Kramers-Kronig is a model-independent transform that checks the EIS data for consistency.
The **Kramers-Kronig** window appears.
4. To clear all fits from the plot, choose **Clear All Fits**.
5. To change time or impedance units, choose **Options**.
This option let you normalize the data and fits to the normalized area.













The Model Editor

The **Impedance Model Editor** allows you to create an equivalent circuit, via a drag-and-drop method.



Don't FORGET! There are several pre-loaded models. Often users find it convenient to start with one of these models and edit it as needed.

Circuit Elements

Symbol	Element	Comments
	Resistor	Abbreviated as R. $Z = R$
	Capacitor	Abbreviated as C. $Z = -i/\omega C$
	Inductor	Abbreviated as L. $Z = iL\omega$
	Constant Phase Element	Models an inhomogeneous property of the system, or a property with a distribution of values. Often abbreviated as CPE.
	Wire	Connects one element to the next.
	Gerischer element	Models a reaction in the surrounding solution that happened already; also used for modeling a porous electrode. Often abbreviated as G.
	Infinite Warburg	Models a linear diffusion to an infinite planar electrode. Often abbreviated as W.
	Bounded Warburg	Models diffusion within a thin layer of electrolyte, such as electrolyte trapped between a flat electrode and a glass microscope slide. Often abbreviated as T.
	Porous Bounded Warburg	Models diffusion through a thin layer of electrolyte, such as electrolyte trapped between an electrode and a permeable membrane covering it. Often abbreviated as O.
	Bisquert Open (BTO)	Transmission line with $Z_b \equiv \infty$ boundary condition. Only pore reactions are considered; base electrode is insulating.
	Bisquert Short (BTS)	Transmission line with $Z_b \equiv 0$ boundary condition. Only pore reactions are considered; electrode is non-insulating.
	Unified (UTL)	Transmission line with a variety of free parameters.

Building an Equivalent Circuit


1. Adding an element

- a. Click on an element symbol.

The element appears in the central window.

- b. Place the mouse cursor over the element. Left-click and drag to move the element to its desired position.

2. Connecting elements

- a. Click on the **Wire** symbol .
- b. Left-click one end of the wire and drag the end to the element.

The element's border turns green when the wire's end reaches the element.

3. Deleting an element


- a. Right-click on the element.

The **Delete** command  appears.

- b. Left-click on the **Delete** command.

The element vanishes.

Here is an example of a simple equivalent circuit (a Randles model) constructed in the **Impedance Model Editor**:



4. Relabeling and fixing parameters for an element



This lets you rename the element, and specify a Lower and Upper Limit for its value. Renaming the element helps you distinguish between elements of the same type during fitting. Giving the program limits on the parameters may help the mathematical algorithm. For example, we know values are generally positive, so a Lower Limit = 0 is reasonable to set.

- a. Left-click on the name of the element (here, R1).

The **Parameter** window appears.



- b. Enter a new **Parameter Name**.
- c. Enter an **Initial Value**, i.e., the first trial value for fitting.
- d. In the **Lower Limit Test** and **Upper Limit Test** fields, enter lower and upper limits, and check the **Enable** checkbox, as desired.
- e. Click the **OK** button.

The **Parameter** window closes, and the element is set to these parameters.


Compiling the Equivalent Circuit


When the equivalent circuit is complete, the circuit can be compiled before use to check for connectivity of the wires. Compiling is only used to check connections

1. Click the **Test Compile** button  in the toolbar. 
The software compiles the equivalent circuit.
If there is a problem, such as a missing connection, an error message appears, and a red box outlines the problem element:
2. Click the **OK** button to continue.
3. Inspect the schematic and make necessary corrections.



If the equivalent circuit compiles properly, the **Model Editor** window appears:


4. Click the **OK** button to continue.
5. Save the equivalent circuit with a *.mdl extension by clicking the **Save** button .



The **File Save As** window appears.


The default folder for saving model equivalent circuits is the **Models** folder.

6. Name and save the file here, or choose a different folder.
The model shown above was saved as `trial model.mdl`.
The **File Save As** window closes.



Fitting the Data to the Equivalent-Circuit Model


1. With the data open and plotted, click **Impedance**, and choose **Fit A Model (Simplex Method)**.



The **Select Model File** window appears.

2. Choose the desired model.

The default folder for models is the Models folder. This Models folder is in the C:\Documents and Settings\All Users\Application Data\Gamry Instruments\Echem Analyst\Models by default. As our example, we choose the model trial model.mdl created previously.




3. Click the **Open** button.

The **Select Model File** window closes, and the **Impedance Fit by the Simplex Method** window appears.

4. Click the **AutoFit** button, or set parameters.

When setting parameters, choose the **Maximum Iterations** to loop before stopping the fit. Enter estimates for all the circuit elements in the **Model Parameters** area. Fix particular elements by enabling their **Lock** checkboxes.



In our example, we try 100 Ω for R_u , 2500 Ω for R_p , and 100 nF for C_f and leave all of them free (unlocked).

- If you don't use the **AutoFit** button, click the **Calculate** button to start the fit. The software attempts to fit the model to the data. When finished, the fitted parameters appear next to each circuit element.


Our model results give

$$R_p = 3 \text{ k}\Omega$$


$$R_{\text{solution}} = 199.7 \text{ }\Omega$$

$$C_f = 980 \text{ nF}$$


Don't FORGET!
If your fit fails, try different starting **Model Parameters**, or even a different equivalent-circuit model.



Like other Echem Analyst fits, the fit also appears superimposed upon the data and a new tab is created that contains those results.



If you try another fit using the same model, this fit will be overwritten. If you fit to another model, the fit results of both models will be displayed.



This new tab shows the residual errors and goodness of fit, along with the various plotting tools. Residuals are a point-by-point **Goodness of Fit**, which quantifies how closely the data match the fit. A smaller number indicates a better fit.

The blue data (**Zreal**) correspond to the y1-axis (on the left); the green data (**Zimag**) correspond to the y2-axis (on the right).

Headings in Data-File Columns

DC Data Files

Abbreviation	Meaning
Pt	Point number
T	Time
Vm, Vf	Measured voltage
Im	Measured current
Vu	Uncompensated voltage
Sig	Signal from the signal generator
Ach	Auxiliary channel
IE Range	I/E (Current Measurement) range on which measurement was made
Over	Any overloads. Numeric record of different overload types
0	No overloads

EIS Data Files

Abbreviation	Meaning
Freq	Frequency
Zreal, Zimag, Zmod, Zphz	Calculated values of impedance
Idc, Vdc	DC component of current and voltage,
Yreal, Yimag	Admittance (calculated from Z)

Current Conventions According to Framework™ and Echem Analyst

The current convention in the Framework software for all experimental packages is that an anodic/oxidation current is positive.

To change the current convention (whether anodic/oxidation currents or cathodic/reduction currents are positive), in the menu **Tools > Options > Units** tab, specify the current you want represented as positive. The current convention can be changed by editing the experimental script (contact Gamry Instruments or your Gamry Instruments representative if you need to do this). Regardless of the current convention used in the Framework, it can be changed in the Echem Analyst to the one you desire (see below for exceptions).



The current convention affects all experiments run under the Physical Electrochemistry and Pulse Voltammetry heading. No other data files are affected.

To change the current convention in the Echem Analyst, in the menu **Tools > Options > Units** tab specify the current you want represented as positive.


To change the current convention in other experimental packages (DC Corrosion, EIS, etc.) please contact Gamry Instruments or your Gamry Instruments representative.

To Edit Visual Basic Scripts:

1. In the toolbar, choose **Tools..**
A drop-down menu appears.



Echem Analyst runs on “Open Source” scripts written in VBA. Most customized analysis routines are done by Gamry Instruments in the factory for you, the user, and that makes Echem Analyst extremely flexible. The typical user will never need to edit the scripts for electrochemical analysis.




Simulating an EIS Curve

It is often useful to simulate the response of an equivalent circuit.


1. Launch the Echem Analyst.
2. Select **Tools > Run Named Script**. Select EIS Model Simulation Script.Gscript.

This opens a blank chart.




3. Select **EIS Simulation > Simulate** (use the Model Editor to build or edit the model).

The **EIS Model Simulation** window appears.




4. Select the saved model, and input parameters for the experiment (frequencies and data-point density) and values of all circuit elements.
5. Click the **Simulate** button.

The simulation appears under new tabs.



This is a simulated Bode plot.



This is a simulated Nyquist plot.

Index

	*		
*.DTA		7	
*.GData		7	
*.mdl		36	
	3		
3D/2D		10	
	A		
AC Calibration Date		22	
Accept		26	
Ach		40	
Ach Range		22	
Ach select		21	
Active Trace area		12	
Add E Constant		17	
Add I Constant		17	
Analysis menu		7	
Apply Template		11	
Apply User-Defined Chart Template window		11	
As a Bitmap		14	
As Text		14	
ASCII		7, 14	
AutoFit button		37	
Automatic Baseline		25, 26	
Automatic radio button		29	
	B		
Beta A		28	
Beta C		28	
bitmap image		14	
Bode tab		32	
Bounded Warburg		34	
	C		
C from CPE, omega(max)		17	
C from CPE, R(parallel)		18	
Cable ID		22	
Calculate button		28, 30, 31, 38	
Capacitor		34	
CD-ROM		7	
Changing the axes		12	
Clear All Fits		33	
Clear Lines		25	
Clear Peaks		25	
Clear Regions		25	
Close button		33	
Color		10	
Common Tools menu		15, 17	
computer		3, 4	
Conditioning		19	
Constant Phase Element		34	
contract engineering		3	
Control Amp Speed		21	
Control Mode		21	
Copy to Clipboard button		10, 14	
Corrosion Rate		28	
Current Convention		22	
Current conventions		41	
Curve Selector		12	
Curve Selector area		11, 12	
Curve Selector button		12	
Cutting and pasting		14	
Cycles		19	
Cyclic Voltammetry		24, 27	
	D		
Data Viewer		10	
DC Calibration Date		21	
Delete		26, 35	
Delta Ep		25	
DigiElch		25	
Draw a Freehand Line button		26	
Draw Freehand Line		11	
	E		
E Log I Fit window		31	
E Log I tab		31	
Ecorr		30, 31	
EIS		32	
EIS Model Simulation window		42	
EIS Simulation		42	
Electrode Area		19	
element		35	
Element		33	
Enable checkbox		35	
Equil. Time		19	
equivalent circuit		33, 34, 35, 36, 42	
Excel®		7, 14	
Experimental Notes tab		8, 20	
Experimental Setup tab		8, 19	
Export to DigiElch		25	
	F		
File Save As window		36	
Final E		19	
Firmware Version		22	
Fit A Model (Levenberg-Marquardt Method)		33	
Fit A Model (Simplex Method)		33, 37	
Framework Version		21	
Framework™ software		6, 24, 41	
Freq		40	
	G		
Gallery		10	
Gamry Analysis Framework Options window		8	
Gamry Echem Analyst Help window		16	
Gamry Help		16	
GamryChart Properties window		10	
General tab		8	
General Toolbar		9, 10	
Gerischer element		34	

Goodness of Fit 39

H

Hardware Settings tab 8, 21, 31
help 16
Help 16
Horizontal Grid 10

I

I/E AutoRange 21
I/E Range 22
I/E Range Lower Limit 21
I/E Range Mode 19
I/E Stability 22
Ich AutoRange 21
Ich Filter 21
Ich Offset 21
Ich Offset Enable 21
Ich Range 21
Idc 40
IE Range 40
Im 40
Impedance Fit by the Simplex Method window 37
Impedance menu 33
Impedance Model Editor 33, 34, 35
Impedance Model Editor window 33
Impedance Subtraction window 33
Inductor 34
Infinite Warburg 34
Init. Delay 19
initial delay 8
Initial Delay 19
Initial E 19
Initial Value 35
Installation 7
Instrument Version 22
Integrate 25, 26, 27
Integrating the voltammogram 26
IR Comp 19

K

Kramers-Kronig 33
Kramers-Kronig method 33
Kramers-Kronig transforms 33
Kramers-Kronig window 33

L

Legend Bar 10
Levenberg-Marquardt method 33
Linear Fit 17, 18
Lock checkbox 37
Log Scale 15
Lower Limit Test 35

M

Mark Found Peaks 11
Max Current 19
Maximum Iterations 37

Min/Max 25, 28
Model 33
Model Editor 33, 34, 42
Model Editor window 36
Model Parameters area 37
Models folder 36, 37
Mouse button 13
My Gamry Data folder 7, 8

N

Normalize by Scan Rate 25
Normalize by Square Root of the Scan Rate 25
Notes... field 20
Nyquist tab 32

O

OK button 33, 35, 36
open 7
Open button 37
Open Circuit 19
Open Circuit Voltage tab 8, 23
Open function 7
Open With... 14
Options 8, 25, 28, 29, 33
Origin® 7, 14
Over 40

P

Page setup window 10
Palette button 10
PaletteBar 10
Parameter Name 35
Parameter window 35
Path 8
Peak Baselines 25
Peak Find 25
Peak Find function 13
plots 9
Polarization Resistance 17, 18, 28, 29, 31
Polarization Resistance window 28
Porous Bounded Warburg 34
Positive Feedback IR Comp 21
Positive Feedback Resistance 22
Post-Run iR Correction 18
Potentiodynamic 30
Potentiostat 21
Powerpoint® 14
Print 10
Print preview 10
Properties... 10
Pstat Model 21
Pt 40

Q

Quick-Integrate 11, 25, 26, 28

R			
Randles model	35		
Region Baseline	27		
Region Baseline Settings window	27		
Region Baselines	25, 27		
Resistor	34		
Reversed	15		
Rotate	10		
Run Named Script	42		
S			
Sampling Mode	19		
Save as Defaults button	29		
Save button	36		
Save Template	11		
Save User-Defined Chart Template window	11		
Scan Rate	19		
Seed Values area	28, 30		
Select Model File window	33, 37		
Select Portion of Curve using the Keyboard	11		
Select Portion of Curve using the Mouse	11		
Select x region	11		
Select y region	11		
Selecting portions of a curve	13		
Selection toolbar	9, 11, 13		
service contract	3		
Show curve selector	11		
Show Disabled Points	11		
Sig	40		
Simplex method	33		
Simulate	42		
Simulate button	42		
Simulating an EIS curve	42		
Smooth Data	17, 18		
Spectrum	33		
Step Size	19		
Subtract Background from File	25		
Subtract Impedance...	33		
support	3		
T			
T	40		
Tafel constants	28, 29		
Tafel Fit	13, 30		
Tafel Fit function	13		
Tafel Fit window	30		
Tafel tab	31		
telephone assistance	3		
Test Compile button	36		
Test Identifier	19		
		Time	19
		toolbars	9
		Tools menu	10, 41
		Transform Axes	18
		Transform Axes selection	15
		U	
		Units tab	41
		Upper Limit Test	35
		Use Seed Values	30
		V	
		Vch AutoRange	22
		Vch Filter	22
		Vch Offset	22
		Vch Offset Enable	22
		Vch Range	22
		Vdc	40
		Vertical Grid	10
		Vf	40
		Visible Traces area	12
		Visual Basic	41
		Vm	40
		Vu	40
		W	
		website	7
		Wire	34, 35
		Word®	14
		X	
		X-Axis column	12
		Y	
		Y2-Axis column	12
		Y-Axis column	12
		Yimag	40
		Yreal	40
		Z	
		Z-clustered	10
		Zimag	39, 40
		Zmod	40
		Zoom	10
		Zphz	40
		Zreal	39, 40