

Redefining Electrochemical Measurement






The Gamry Instruments Flexcell[™] is a versatile electrochemical cell designed to counteract a problem that plagues most cells -- crevice corrosion around the specimen seal. First popularized by Avesta Steel, the cell utilizes a flooded gasket seal design to inhibit crevice corrosion between the sample and its holder. This simple yet elegant design results in an easy to use, reliable, crevice-free system.

The main body of the cell is a two piece construction that allows for easy handling and assembly. Specimens are mounted facing upwards, allowing for close but unobtrusive placement of a temperature probe and a reference electrode. The Flexcell can be heated externally with a Heating Mantle, internally with insertion heaters, or with a Recirculating Bath. A glass or Teflon-coated copper coil is supplied for temperature control of the cell using a Recirculating Bath.

When used in conjunction with CPTIIO[™] Critical Pitting Temperature Software and a TDC2[™] Temperature Controller, the Flexcell becomes part of a totally automated tool for the determination of Critical Pitting Temperature (ASTM G I50). In the tradition of Gamry, the Flexcell can be used out of the box or adapted to meet your needs.



Cross Section View of the Flooded Gasket Seal

Specifications

| Unless otherwise specified, these specifications refer to the flushed port configuration. | |
|---|--|
| Vessel | |
| Volume | 1000 ml (minimum) 1300 ml (maximum) |
| Inner diameter | 10.0 cm (nominal) |
| Metal Specimen | |
| Active Area | 0.775 square inches (nominal) 5.0 cm ² (nominal) |
| Minimum Size | 4.0 cm x 4.0 cm 4.0 cm diameter (circular) |
| Maximum Size | 7.5 cm x any length or 11.0 cm x 11.0 cm (rectangular) 11.0 cm diameter (circular) |
| Temperature Range | |
| Operating | -25 °C to 140 °C |
| Storage | -40 °C to 180 °C (not assembled) |
| Chemical Compatibility | |
| Wetted Materials | Pyrex, Teflon, Filter Paper |
| Non-Wetted Materials | Above, plus stainless steel, Viton, silicone rubber |
| | |

Parts Supplied

The CPT Flexcell is shipped with all the accessories and manuals required to run crevice free tests except for these items which are subject to unique customer requirements:

- No metal samples are supplied.
- A reference electrode is not included.

Accessories

For ASTM G I50, "Electrochemical Critical Pitting Temperature Testing of Stainless Steels", the Flexcell should be used with the TDC2 Temperature Controller and the Heating Mantle (Part No. 990-I47/8). For convenient and effective stirring of the Flexcell, use the Flexcell Stirring Rod Kit (Part No. 990-II7) and Stirring Motor Kit (Part No. 990-I52/I53).

Rev 6.0 8/19/2011 ©Copyright 1990-2011 Gamry Instruments, Inc. All specifications subject to change without notice.



+215-682-9330

Warminster, PA 18974

Fax: +215-682-9331