



Load/Power Supply Interface (LPI1010) High Voltage EIS Test System QUICK-START GUIDE



The full experimental setup requires the following:

- Gamry Instruments Interface™ 1010 Potentiostat/ Galvanostat/ZRA
- An LPI 1010 Load/Programmable Power Supply Interface containing:
 - An LPI Top Assembly
 - An LPI Power Cable
 - An LPI Voltage Sense Cable Kit
 - Gamry Instruments Framework™ Software

Install your IFC1010 according to the IFC1010Interface 1010 UserOperation Manual and Quick -Start Guide available at :

<https://www.gamry.com/support/documentation-downloads/>.

Ensure that your Device Under Test is not connected to your setup.

1

Connect the LPI Top Assembly to the Interface 1010 and LPI Power Cable.

2




Plug LPI Top Assembly into the Cell Cable receptacle.



Plug LPI Power Cable into the small receptacle on the top of the LPI Top Assembly.

3

Connect the LPI Power Cable to the User I/O Connector on the rear panel of the Interface 1010.




Rear panel of Interface 1010

LPI Power Cable

SUB Module



- ④ Connect the Device Under Test (DUT) and the LPI Top Assembly to your Electronic Load or Power Supply.




5



Connect the DUT and the Cable End Modules together via the LPI Voltage Sense Cables.



- ⑥ Configure your Electronic Load or Programmable Power Supply for **current control** and set the full-scale current as appropriate for your experiment.
- ⑦ Launch Framework™ software and open the Gamry Instrument Manager software.
 - Select the appropriate instrument on the left hand side of the instrument manager.
 - Press the Configure button as shown.



- Select a power device from the list of available power devices. The rest of the fields should autopopulate.



If your device is not in the list, you need a Power Device configuration file. Please contact your Gamry representative to discuss your options.

WHAT DOES GAMRY SOFTWARE DO?



Gamry Framework™

Potentiostat control for flexible data acquisition. Select from standardized experiments grouped by research type, or use the Sequence Wizard to build complex automated experiments.



Echem Analyst™

Quick and easy data analysis. Open data files with Echem Analyst for specialized analysis algorithms and high-quality plots. Customize, overlay, and scale plots, or export data.



My Gamry Data™

The default data-folder location for Gamry Framework, with a shortcut on your desktop after installation. Change the folder location within Gamry Framework via **Options > Path**.



Virtual Front Panel™

Software-based front panel for quick access to Gamry potentiostats' functions, like a front panel of an early analog potentiostat; and to perform simple electrochemical experiments.



Electrochemical Signal Analyzer™

Designed specifically for the acquisition and analysis of time-dependent electrochemical noise signals.



Resonator™

Data-acquisition and -control software for the Gamry eQCM™. Contains a full suite of physical electrochemistry techniques.



Electrochemistry Toolkit™

A sophisticated package for complete access to the capabilities of Gamry potentiostats in the software environment of your choice.