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*Evaluation of Organic Coatings with Electrochemical Impedance Spectroscopy.*

*Part 1: Fundamentals of Electrochemical Impedance Spectroscopy.* David Loveday, Pete Peterson, and Bob Rodgers, *JCT CoatingsTech*, 46-52, August 2004.

*Evaluation of Organic Coatings with Electrochemical Impedance Spectroscopy.*

*Part 2: Application of EIS to Coatings..* David Loveday, Pete Peterson, and Bob Rodgers, *JCT CoatingsTech*, 88-93, October 2004.

*Evaluation of Organic Coatings with Electrochemical Impedance Spectroscopy.*

*Part 3: Protocols for Testing Coatings with EIS.* David Loveday, Pete Peterson, and Bob Rodgers, *JCT CoatingsTech*, 22-27, February 2005.

*Electrochemical Impedance and Noise,* Robert Cottis and Stephen Turgoose, NACE International, 1440 South Creek Drive, Houston, TX 77084, US, 1999.

Website: [www.nace.org](http://www.nace.org).

*Impedance Spectroscopy; Theory, Experiment, and Applications,* 2nd ed. , E. Barsoukov, J.R. Macdonald, eds., Wiley Interscience Publications, 2005.

*Electrochemical Methods; Fundamentals and Applications,* A.J. Bard, L.R. Faulkner, Wiley Interscience Publications 2000.

*Electrochemical Impedance: Analysis and Interpretation,* J.R. Scully, D.C. Silverman, and M.W. Kendig, editors, ASTM, 1993.

*Physical Chemistry,* P.W. Atkins, Oxford University Press 1990.

*Signals and Systems, A.V. Oppenheim and A.S. Willsky, Prentice-Hall, 1983.*

*The Use of Impedance Measurements in Corrosion Research; The Corrosion Behaviour of Chromium and Iron Chromium Alloys, J.A.L. Dobbelaar, Ph-D thesis TU-Delft 1990.*

*Characterization of Organic Coatings with Impedance Measurements; A study of Coating Structure, Adhesion and Underfilm Corrosion, F. Geenen, Ph-D thesis, TU-Delft 1990.*

*Identification of Electrochemical Processes by Frequency Response Analysis, C. Gabrielle, Solartron Instrumentation Group 1980.*

*Comprehensive Treatise of Electrochemistry; Volume 9 Electrodics: Experimental Techniques; E. Yeager, J.O'M. Bockris, B.E. Conway, S. Sarangapani, chapter 4 "AC Techniques", M. Sluyters-Rehbach, J.H. Sluyters, Plenum Press 1984.*

*Mansfeld, F., "Electrochemical Impedance Spectroscopy (EIS) as a new tool for investigation methods of corrosion protection", Electrochimica Acta, 35 (1990) 1533.*

*Walter, G.W., "A review of impedance plot methods used for corrosion performance analysis of painted metals", Corrosion Science, 26 (1986) 681.*

*Kendig, M., J. Scully, "Basic aspects of electrochemical impedance application for the life prediction of organic coatings on metals", Corrosion, 46 (1990) 22.*

*Fletcher, S., "Tables of Degenerate Electrical Networks for Use in the Equivalent-Circuit Analysis of Electrochemical Systems", J. Electrochem. Soc., 141 (1994) 1823.*

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### Technical Support

A compilation of technical support information by hardware and software.

### Software Updates

Licensed users can download the latest version of the Gamry Software.

### Documentation Downloads

Download manuals for hardware and software products.

### Electrochemistry Courses and Training

Information on short courses and electrochemistry training.

### Contact Technical Support

Get support from the electrochemical experts at Gamry.

### Support & Tutorial Videos

For all of our tutorial and support-related videos you can also visit Gamry Instruments YouTube

### Application Notes

A series of Application and Technical Notes to assist you in getting the most accurate results.

### Literature Database

A database of articles where people have used a Gamry System in their research

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English ▾

# *Market Leader in the Support of Electrochemical Research*

## *Sign Up For Our eChem Newsletter*

Stay up-to-date on products, conferences and news.

NAME

EMAIL ADDRESS

X

S

olled off the  
rical isolation in the


industry. The ability to couple with loads and power supplies is a good indicator of our floating capability.

[LEARN MORE](#)

## *Open-Source Scripting* for easy customization, supporting unique experiments and set ups.

Gamry Instruments has long given you the ability to modify your standard experiments as needed, or even develop your own experiments. All of our experimental scripts are written using an open-source scripting language, developed by Gamry, called Explain™.

[LEARN MORE](#)

## *Higher Accuracy as demonstrated by accuracy contour plots.*

The Accuracy Contour Plot gives you valuable information about a potentiostat's EIS performance and accuracy. It allows you to compare the Electrochemical Impedance Spectroscopy capability of any potentiostat.

[LEARN MORE](#)

## *Deep Expertise Scientists in all areas of operations to ensure*

# quality and efficacy

Our sales and support staff consists of electrochemists and corrosion scientists available to discuss your applications, assist you with data analysis, or just offer advice on how to make the most out of your Gamry Instrument.

[LEARN MORE](#)

## Technical Resources

### APPLICATION NOTES

### PRODUCT SUPPORT

### LITERATURE DATABASE

Gamry has been at the forefront of electrochemical corrosion testing since our inception in 1989. Our customers have told us that ease of use and reduced acquisition time are the most important features for their research.

## Market Applications

Our instruments are suited for battery testing of materials to coin cells and all the way to large format cells. We even have a system for stack testing. All of our potentiostats are electrically isolated from

our paints and coatings customers have told us that ease of use and reduced acquisition time are the most important features for their research.

Photovoltaic researchers looking for a complete solution will find Gamry's IMPS/IMVS quite capable of helping the researcher find the answers they need.

earth-ground.

## CORROSION


Most sensors are amperometric or impedimetric. Gamry makes it easy for researchers to generate and analyze data from either of these sensors.

## SENSORS

## FUEL CELLS

## COATINGS

Testing for supercapacitors is very similar to testing of batteries, just on a faster time scale. We have a number of recommended setups for the supercapacitor researcher.

## SUPERCAPACITORS

## PHOTOVOLTAICS

Our software toolkits let OEMs take control of our instruments from their own software providing a seamless interface for data acquisition.

## OEM / PARTNERS

# Featured Product

## The New Interface 1010 Series

The Interface 1000 potentiostats have been replaced with the new Interface 1010 Series now with built-in temperature monitors using an RTD probe. The Interface 1010E has higher EIS frequencies with extended EIS to 2 MHz.

[GET THE SPECS & ACP](#)

#### HOME OFFICE

734 Louis Drive

Warminster, PA 18974

P: 215-682-9330

---

#### ABOUT US

[Life & Times](#)

[Our Team](#)

[Privacy Policy](#)

[Gamry Difference](#)

#### ONLINE STORE

[Create Account](#)

[Login](#)

[Accessory Quote](#)

[Customer Service](#)

#### APPLICATION NOTES

[Impedance Spectroscopy](#)

[Battery Research](#)

[Corrosion/Coatings](#)

[Research Electrochemistry](#)

#### LATEST NEWS

[General News](#)

[Conferences](#)

[Echem Newsletter Archive](#)

[Careers](#)

#### CONTACT US

[Sales](#)

[Tech Support](#)

[Find a Distributor](#)

[Chat with Support](#)

#### APPLICATIONS

[Battery Testing](#)

[Corrosion](#)

[Fuel Cells](#)

[Paints & Coatings](#)